

## Instruction for use

# Sodium hyaluronate 1 mg/ml HYLO COMOD

To improve eye surface moisture in eyes that are dry, burning or irritated due to environmental conditions.  
HYLO COMOD is a preservative free and phosphate free solution.  
Lubricating eye drops.

### What is HYLO COMOD?

#### Composition:

1 ml of HYLO COMOD contains

Active ingredient: Sodium hyaluronate 1 mg

Other ingredients: Citric acid anhydrous, sodium citrate, sorbitol and water for injections.

### What is HYLO COMOD used for?

To improve the eye surface moisture in eyes that are dry, burning or have foreign body sensations due to environmental conditions.


### When should HYLO COMOD be used?

If your eyes are itching, burning, or easily tiring, if you feel, as if you had sand in your eyes, this might be due to a lack of moisture on the surface of your eyes. If you spent a longer time in an airconditioned room, for example, or after intense computer work the natural tear film that usually covers the eye surface breaks up more quickly. For your eye tissue to be fully functioning, it needs to be continuously lubricated by a tear film.


HYLO COMOD contains the sodium salt of the hyaluronic acid, a natural substance that is found in the eye but also in other parts of the body. Its special property is that it forms a regular, stable, long lasting lubricating moisture film on the eye surface that is not easily rinsed off. It does, however, not cause blurred vision and protects your eyes from dryness and irritation for a long time.

This is particularly important if you wear contact lenses, because the contact lens needs to swim in sufficient tear fluid in order not to cause any discomfort. HYLO COMOD may be used while wearing contact lenses. Wearing soft or hard contact lenses will be made more comfortable by using HYLO COMOD as it does not form crusts or residue.


### How should HYLO COMOD be applied?



Remove cap before use.



Before the first application please turn the HYLO COMOD bottle with its tip down and press onto the bottom of the bottle until the first drop appears at the nozzle. Now the bottle is ready for further use. Hold the bottle upside down and apply a strong and quick pressure on the bottom of the bottle. This activates the mechanism for the release of one drop. Because of the special COMOD valve system, size and speed of each drop are the same however much pressure is applied.



Lean your head back slightly, pull down your lower eyelid and instil one drop into the lower eye. Slowly close your eye to allow the fluid to spread evenly on the eye surface.



Replace cap immediately after using. Please take care that the dropper tip is dry. When applying the fluid please avoid any contact between the tip of the bottle and your eye or your skin.

The COMOD - system releases 10 ml of solution which is the equivalent of approximately 300 drops. For technical reasons a small amount remains in the bottle at the end of exhaustion.

HYLO COMOD can be used for an unlimited period of time.

Only one person should be treated with one particular bottle HYLO COMOD.

### **What else do you need to know?**

If you use any other eye drops there should be at least 30 minutes between the applications and HYLO COMOD should always be used last. Eye ointments should, however, always be administered after using HYLO COMOD.

Because it is preservative-free HYLO COMOD is well tolerated even when used over a longer period of time. In very rare cases hypersensitive reactions like burning or excessive tearing have been reported which stopped immediately when the use of HYLO COMOD was discontinued.

### **How often and for how long can you use HYLO COMOD?**

HYLO COMOD should be used according to the instructions given by your doctor or your contact lens specialist.

### **Posology and administration**

Generally you instil one drop of HYLO COMOD three times a day into each eye. If necessary it can also be used more frequently. If you use HYLO COMOD more frequently (e. g. more than 10 times per day) please consult your ophthalmologist. HYLO COMOD is suitable for longterm treatment.

### **Contraindications**

HYLO COMOD should not be used if you are oversensitive to one of the ingredients.

### **Special warnings and precautions for use**

Do not touch the nozzle and do not allow the nozzle tip to touch the eye during use. Do not allow your HYLO COMOD to be used by anyone else.

Do not apply HYLO COMOD solution at the same time as other eye medications; eye drops can be applied 30 minutes after using HYLO COMOD. Ointments should always be applied after using HYLO COMOD.

### **Interaction with other medicinal products and other forms of interactions**

None.

### **Pregnancy and lactation**

No well-controlled studies in pregnant or breast-feeding women exist concerning sodium hyaluronate. Nevertheless, no negative effects on the fetus or the infant are to be expected.

### **Driving and using machines**

The eye drops may cause blurred vision for a short time after application even at normal dosages and with proper use. This can subsequently impair reaction time while driving or operating machinery.

### **Undesirable effects**

Because it is preservative-free, HYLO COMOD is well tolerated even when used over a longer period of time. In very rare cases, hypersensitive reactions like burning or excessive tearing have been reported which stopped immediately after the use of HYLO COMOD was discontinued.

**Warning:** If irritation persists or increases, discontinue the use and consult the physician.

### **Overdosage**

Given the nature of product and the route of administration, no problems of overdosage are expected as the excess fluid will flow from the eye.

### **Incompatibilities**

Do not use concomitantly with disinfectants containing quaternary ammonium salts because hyaluronic acid can precipitate in their presence.

### **What you need to know about storing HYLO COMOD?**

Shelf life as packaged for sale: 3 years

Store in a cool and dry place.

Keep out of the reach of children.

HYLO COMOD can be used for 12 weeks after opening.

Do not use beyond expiry date.

### **Nature and contents of the container:**

10 ml of eye drops solution in multidose plastic container with airless pump and closed with a cap.

### **Date of information:**

October 2018

### **Imported and marketed in India by:**

URSAPHARM India Pvt. Ltd.

PLOT NO A-1, FIEE, BASEMENT, OKHLA INDL. AREA PHASE II, NEW DELHI-110020

Import License no.: IL/FF-000279 (FF-421)

### **Manufactured by:**

URSAPHARM Arzneimittel GmbH, Industriestraße 35, 66129 Saarbrücken, Germany

Tel.: +49 (6805)/9292-0, Fax: +49 (6805)/9292-88,

email: info@ursapharm.de,

www.ursapharm.de